

**HESTION**

**- Histology**

**Daintree Scientific – Australia**



**LST94 SLIDE STAINER**



**Membrane keypad**

**Microcomputer technology advanced intelligent control system**

**Two slide racks and two programs running simultaneously greatly improve working efficiency. Besides the stainer function, it also has a flexible program adaptable for different processing and lab tests. This flexible program can begin or end at any station as desired , and is capable of up to 25 steps**

**The system can store four different programs which apply to various tissue processing protocols. It also has wash function**

**The system has full protection . Even if any abnormal conditions arise, or if operated incorrectly, it cannot damage the equipment**

**Automatic position search function means the basket can move to the set position from any other position**

**Alarm function monitors the operation being carried out and prevents any incorrect operation in time for the operator correct it.**

**It uses high-performance stepping motor systems to ensure running stability, low noise and accurate positioning**

**Air filter device to reduce air pollution**

**Sealed soft-touch control panel for reliable operation ,dustproof and waterproof**

**The streamlined exterior shape makes it easier to clean**

**It has inlet and outlet water piping to wash the tissue at any time**

**Technical Parameters:**

**Quantity of slide racks :2**

**Throughput capacity: 96 slides at a time(48 slides on each slide rack).**

**Number of containers: 19**

**Volume of container: 780ml for each container, the first and second containers comprise the washing station with circulated water**

**Retention period in each station: 1 second-90 minutes 59 seconds, can be set as desired**

**Environmental protection: suction and filtration of fumes using an activated carbon system which prevents any toxic gases from entering the laboratory environment**

**LCD display: with user-friendly keyboard to program and display all information concerning staining**

**Protection function: it can prevent damage to tissue samples should the basket movement be physically blocked while in operation**

**HESTION**

**Daintree Scientific Tel: (03) 6376 3335 Fax: (03) 6376 3334**

**Address: 20 Kiama Parade, St. Helens, Tas 7216 ,AUSTRALIA**

**Visit our website: [www.daintreescientific.com.au](http://www.daintreescientific.com.au)**

**EMAIL:**

Asia & Middle East :	info@histo-instrument.com
Europe :	hestion-eu@histo-instrument.com
Africa :	hestion-af@histo-instrument.com
USA :	hestion-us@histo-instrument.com
Australia & New Zealand :	info@daintreescientific.com.au