

HESTION

- Histology

Daintree Scientific – Australia



TEC2900 HISTO TRIM

The TEC-2900 Histo Trim is a practical device designed to remove excess paraffin from the embedding cassettes after inclusion of specimens. By scraping the embedded block on the “grooved plane” heated surface, excess paraffin will be removed leaving the cassette clean , ready to be set properly on the clamp of the microtome.

Temperature can be set from ambient to 90°C

Table top (mm): 150×220

Temperature control accuracy: ±1°C

Rated voltage: 220V ± 10% 50HZ

Outside dimensions (mm): 200(w)×160(D)×175(H)

Weight: 2.5kg

HESTION

Daintree Scientific Tel: (03) 6376 3335 Fax: (03) 6376 3334

Address: 20 Kiama Parade, St. Helens, Tas 7216 ,AUSTRALIA

Visit our website: www.daintreescientific.com.au

EMAIL:	Asia & Middle East :	info@histo-instrument.com
	Europe :	hestion-eu@histo-instrument.com
	Africa :	hestion-af@histo-instrument.com
	USA :	hestion-us@histo-instrument.com
	Australia & New Zealand :	info@daintreescientific.com.au